

CS Series Centrifuges

CS1385 / CS2385 / CS1465 / CS2465 — Precision. Performance. Reliability.

The CS Series centrifuges deliver exceptional performance in a compact, reliable package. Whether you need a fixed-speed workhorse or a fully programmable variable-speed system, the CS Series has a model to match your laboratory requirements. All models share the same compact footprint and universal power supply, making them ideal for busy benchtops worldwide.



CS1385 / CS2385

Fixed Angle Rotor | Fixed & Variable Speed

- 8 - Place Fixed Angle Rotor
- 1,712 × g at 3,400 rpm (CS1385)
- Variable speed: 500–4,000 rpm
1,564 × g at 3,400 rpm and 2,164 × g at 4,000 rpm (CS2385)
- Speed adjustable during operation
- Compatible with $\Phi 16 \times 125/115/110/100$ mm tubes (direct fit) and $\Phi 13 \times 75$ mm / $\Phi 10 \times 40$ mm (with adapter)
- Segment LCD display with key-operated panel

CS1465 / CS2465

Horizontal Rotor | Fixed & Variable Speed

- 6 - Place Horizontal Rotor
- Fixed speed at 3,400 rpm with 1,712 × g RCF (CS1465)
- Variable speed: 500–4,000 rpm
2,370 × g at 4,000 rpm (CS2465)
- Compatible with $\Phi 16 \times 110/100$ mm tubes (direct fit) and $\Phi 13 \times 75$ mm / $\Phi 10 \times 40$ mm (with adapter)
- Segment LCD display with key-operated panel

Key Features

High Performance

Up to 4,000 rpm &
2,164 g RCF

Versatile

Fixed Angle and
Horizontal Rotor Options

Safe & Silent

Lid-lock safety &
quiet operation

Intuitive

LCD & touchscreen
models available